



Advance Life Support Ambulance



AMBULANCE



SALIENT FEATURES



- Modular Body
- Plush ABS/FRP Formed Interiors
- Duct Air Conditioner
- Liberal Head Room
- Ample Storage Space
- Ergonomic Seating
- Retractable Doctor's Seat
- Centralised Suction
- Adequate Lighting

We at **ASH** Helathcare design ambulance which is capable to support advance life support to patients allowing transportation at long distances giving them a state of the art facility which is almost identical to the Hospital ICU.


The vehicle interior is designed with a modular effect keeping in mind the provisions required for setting up the equipments that would be required in a Criticalcare environment. All provisions are designed & made available in such a way that it allows user to equip the vehicle with specialised life support equipments like Transport Ventilator, Patient Monitor, Defibrillator, Syringe & Infusion Pump & suction Apparatus. Each equipment provision & placement is done keeping in mind the user convenience for use in the vehicle. Snap Fit arrangements are provided to firmly hold the equipments at the same time allowing easy removal from the vehicle at the time of emergencies or while patient transfer.

The standard Ambulance supplied comes equipped with Central Gas Supply System which has an ON/OFF system & electronic Alarm System in Case of Low Gas. Features like Wash Basin with Auto Water Flow, Waste Bins, Drawers, & storage facilities are offered as standard part of the system & are conveniently placed for easy use.

Apart from the above as standard we offer imported lightweight Automatic Loading Patient Stretcher which can be handled by single person, All Aluminium Foldable Sccop Stretcher & high Intensity imported light Bar with Hooter & public address System.



MULTIPARAMETER PATIENT MONITOR



EVA MP 12

- 12.1" TFT Display offering Continuous monitoring of vital parameters like ECG+SpO2+NIBP+Resp+Temp.
- Suitable for Transport with built-in battery backup


TRANSPORT VENTILATOR



Flight 50

- Compact and Light weight,
- Suitable for Adult & Pediatric
- Invasive & Non Invasive Facility
- Ventilatory modes: A/CMV, VIM, VPIBT (AP, BiPAP)
- Pressure Support, Built in PEEP
- Backup ventilation in all modes

PULSE OXIMETER



IP 1010

- Large LCD display for SpO2, Pulse rate & Plethysmographic waveform
- Easy to operate with AC & Battery Backup
- Audio & Visual alarms.

TRANSPORT CUM ICU VENTILATOR



Flight 60

- Light weight & portable
- Suitable for Adult & Pediatric Use
- High Performance Ventilator with inbuilt Piston
- Ventilation Modes: CPAP, Bilevel S, Bi Level ST, PCV, Asst PCV, PSV, VT+PSV, SIMV, VCV.
- Inbuilt Battery Back-up 12 hours

SYRINGE INFUSION PUMP



Graseby 2000 & 2001

- Compact, Light weight & Portable with inbuilt rechargeable battery
- Suitable for infusion of emergency drugs
- Accepts all the major brands of Syringes.

RESUSCITATION KIT



- Light weight, durable & simple Ventilator
- Ventilation Modes : IPPV, A/C, Manual & Peep
- Kit Includes Ventilator, O2 Cylinder, regulator, face Mask, Laryngoscope, Suction attachment & Carry case

SUCTION MACHINE



MC 600

- Compact, Portable & Light weight
- Inbuilt Rechargeable Battery operated
- High efficiency of motor suction pump.

DEFIBRILLATOR



CU ER5

- Biphasic defibrillator,
- AED Mode with Voice Command, Battery Backup, ECG Monitor, 5 Lead ECG Cable, Light Weight & Compact.



its state of the art Manufacturing unit is growing as a major player providing Healthcare on wheels.

ASH healthcare being a major supplier of Criticalcare equipments now is able to provide an economic solution & a one stop facility where it allows Hospitals & healthcare providers to just look at one supplier who can provide them high quality equipments with can be seamlessly integrated in a vehicle.

We not only design & supplies high quality Ambulance systems but also has designed & supplied other on the wheels healthcare solutions like Mobile Diagnostic Centers, Mobile Dental Clinics, Mobile Operation theaters, Mobile Lithotripsy center & so on.

We at ASH healthcare are poised to offer as many mobile healthcare solutions as possible & we design, develop Vehicles as per user needs & requirements. Whatever you needs be we can always provide you an effective mobile solution enabling you to expand your reach & services to the community at large.



Mobile Dental Clinic



Marketed & Service by :



+91 9925063777 | +91 79 65401100
 sales.india@aussin.co
 www.aussin.co | www.ashbme.com

Authorised Representative

