


OXYCORR

Portable Pulse Oximeter



Features

- OXYCORR provides SpO₂ and pulse rate measurements of any patient ranging from Adult to Neonates
- Option to display Plethysmograph or values in extra big Digits.
- SAC Technology : Enhanced performance in low Perfusion & Artifacts.
- Audible pitch change with change in % SpO₂
- Non-volatile memory for storing Graphical & Tabular trends with alarm status for 72 Hours.
- Serial interface for connecting printer / PC for printing SpO₂ & PR values every 1 Min, 5 Min, 15 Min, 30 min
- Compatible to connect Complete Care CMS & PC-S softwares.
- Long battery operation using power saving mode.
- Display Size: 5.6"